Ally® Wafer Extenders

en Instruction for use

1. Product description: Elastic tape is designed to prevent the edges of the baseplate from lifting and rolling and keeps your barrier securely in place. For use with colostomies, ileostomies or urostomies.

2. Composition: Release Liner, hydrocolloid (CMC, medical hot-melt glue) and backing material.

3. Specification: R80mm

4. Intended purpose: Elastic tape is designed to stick around the edges of your baseplate, to prevent the edges of the baseplate from lifting and rolling and keeps your baseplate securely in place. It is elastic so it follows body contours and movements. It can absorb moisture, and it is also skin-friendly.

5. Intended population: The product may be used on individuals with a colostomy, ileostomy or urostomy.

6. Intended user: Health care professionals, caregivers and patients in a hospital or home environment.

7. Indications: For use with a colostomy, ileostomy or urostomy.

8. Warnings:

8.1 Single use.

8.2 For external use only.

8.3 Do not apply to wounds or damaged skin.

8.4 Should skin irritation occur, contact your doctor or another health professional.

8.5 Do not use on individuals with a known sensitivity to the ingredients.

8.6 Do not use on individuals with impaired immune systems or people receiving immune suppressive treatment.

8.7 Do not use on premature babies, or infants receiving incubator care.

8.8 Do not use if the product expires.

9. Information for use:

9.1 Preparation:

9.1.1 Begin application by warming the elastic tape briefly between your hands to aid adhesion.

9.1.2 Ensure that the skin around the stoma is clean and dry.

9.2 Application:

9.2.1 Remove half of the backing film and apply the elastic tape to the skin, overlapping the edge of the barrier. The elastic tape should overlap the edges of the skin barrier at least half an inch.

9.2.2 When application, do not stretch the product.

9.2.3 Remove the remaining backing and smooth down to aid in adhesion.

10. Method of removal: The elastic tape may come off together with the baseplate. If needed, loosen the edges of the tape and gently stretch the tape for removal. An adhesive remover may be used to assist.

11. Storage: Keep away from sunlight. Store products at cool and dry conditions – do not freeze.

12. Disposal: The product should be disposed of together with household waste. Do not flush the product down the toilet.

13. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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