



2-Piece System

FLEXIMA® KEY

THE KEY TO FEELING SECURE WITH EXTRA FLEXIBILITY

Flexima® Key The guided positioning A major step towards security!

Key Guided Positioning

- O Precise and simple positioning with a guiding system
- Secure set-up decreases risks of leakages
- O Minimum patient education: intuitive coupling



Key guiding system

Secured Coupling System

- O Large stable track, no distorsion.
- Specific adhesive foam formula ensures optimum leak proofing



Pouch fitting

Easy Pouch Removal

- Special cut-away section at the top of the flange
- O Pouch ring leaves no residue when the pouch is removed



Removal

Extra flexibility

Soft material for flexibility

Follows body contours

Smooth and low profile track

• For an easy cleaning and great discretion

Soft convexity to reduce leakage

- For flush or slightly retracted stomas
- A belt can be used for more security



Flexima® Key Drainable Roll'Up pouches


High performance filter

- Water and oil resistant filter No risk of staining clothes Low risk of pouch ballooning





480 ml**



NEW Shape and features

• Ergonomic shape for more wearing comfort Optimized sizing for optimized capacity

Integrated soft Roll'Up clamp

Intuitive and secure closing system




O Easy to use outlet

Mini:

250ml*

O Soft edges for an increased wearing comfort once the Roll'Up is closed












^{**}for diameter 80mm

Flexima® Key Closed and Drainable pouches

Closed pouch



• Laminar technology filter:

 High performances of degassing and deodorizing

Maxi drainable pouch



1. Comfort:

- Funnel-shaped for a perfect output control
- Large diameter for easy emptying
 & direct connection to a Flow Collector
 bag

2. Security:

- Easy to open and close
- Suitable for patients with low dexterity

LKP Filter Deodorized gases Liquid Gases

Water and oil resistant filter

Flexima[®] Key Right positioning – Fit from the start





Very easy to use for a lot of people»
Jos- ET nurse, Netherlands

The Key is a good idea because the problem we're facing with this kind of product is being able to position it correctly.

With Key, you can't go wrong!»

Isabelle- ET nurse, Netherlands

Flexima® Key

Closed pouches (box of 30)	Diameter (mm)	References beige	References transparent
NEW	40	62142	O 62242
	40	62140	O 62240
	50	62150	O 62250
	60	62160	O 62260
	80	62180	○ 62280

Flat base plates (box of 5)		Diameter (mm)	References
Cut-to-fit	×	40/15-35	62040
	*	50/15-45	62050
	×	60/15-55	62060
	*	80/15-75	62080
Pre-cut		40/20	620421
		40/25	620425
		50/30	620430
		50/35	620436
		50/30	620430

Convex base plates (box of 5)		Diameter (mm)	References
Cut-to-fit	*	50/15-25	62052
	×	60/15-35	62062
	*	60/15-45	62063
Pre-cut		50/20	62054
		60/25	62064
		60/30	62066
		60/35	62067
		60/40	62069
Tre-cut		60/25 60/30 60/35	62064 62066 62067

Roll'Up pouches (box of 30)	Diameter (mm)	References beige	References transparent
Mini	40	621243	O 621242
Midi	40	6 21241	O 621240
	50	6 21251	O 621250
Maxi	60	6 21261	O 621260
	80	621281	O 621280

Maxi pouches (box of 30)	Diameter (mm)	References beige	References transparent
	50	• 62550	o 62650
	60	• 62560	o 62660
AVAILABLE SOON	80	62580	o 62680

Stomacare belt	References	
1 unit	5028	

This document, its contents, including institutional data, information, trademarks and logos mentioned herein are the exclusive property of B. Braun. Any representation and / or reproduction, total or partial, of this document and its contents without the express prior consent of B. Braun, is strictly prohibited and constitutes an infringement of the intellectual property rights of B. Braun.

Commercialized and manufactured by:

B. BRAUN MEDICAL SAS | 26 Rue Armengaud | 92210 Saint-Cloud - France

Tel. 01 41 10 53 00 | Fax 01 41 10 75 68 | www.bbraun.com

B. Braun Medical SAS with a share capital of 31 000 000 € | RCS Nanterre 562050856



For most countries a specific letter is added at the end of the radical reference.

All volumes are measured according to ISO 8670-2 standard.

This medical device is a regulated health care product and on the basis of this regulation is certified CE.

The products commercialized by B. BRAUN MEDICAL shall exclusively be used in accordance with the package inserts which have been updated and included in the boxes. Non binding documents and photographs